

Cryotainer B.V.  
 Fokkerstraat 507  
 3125 BD Schiedam  
 Phone: +31 (0)10 262 30 53

P.O. Box 130  
 3100 AC Schiedam  
 www.cryotainer.com



## Datasheet FT 1700 (Flying Tainer)



Dimensions l x h x b	mm	12.200 x 3.650 x 2.500
Maximum working pressure	bar	10,0 or 22,0 ( depending on storage container )
Maximum flow	Nm <sup>3</sup> /h	1.700 ( 2 x 850 ) for 8 hours per day maximum
Maximum outlet temperature nitrogen	°C	15 <sup>0</sup> till 20 <sup>0</sup> below ambient temperature
Total weight	kg	38.000
Suitable for cryogenic gasses		LIN; LOX*
Storage capacity	ltr	19.700
Outlet connection		Both sides 2" weco fig. 100

\* For use of LOX container and vaporizer should be purged before use